

DNA precipitation (To concentrate DNA sample)

- Reagents:
 - 3 M sodium acetate pH 5.2 or 5 M Ammonium or Sodium acetate
 - DNA solution
 - 100% ethanol (ice cold) or -20 degree Celsius
 - 70% ethanol kept on 4 degree Celsius
- Procedure:
 1. Measure the volume of the DNA sample.
 2. Add 1/10 volume of sodium or potassium acetate, pH 5.2
 - These amounts assume that the DNA is in TE only; if DNA is in a solution containing salt, adjust salt accordingly to achieve the correct final concentration.
 3. Mix well vortex.
 4. Add 2 volumes of ice cold 100% ethanol.
 5. Mix well by vortex.
 6. Place on ice or at -80 degrees C for 30 minutes. Or at -20 overnight.
 7. Spin a maximum speed in a microfuge 10-15 min.
Check for DNA pellet little opaque type
 8. Carefully decant supernatant.
 9. Add 500ul 70% ethanol. Mix gently. Carefully decant supernatant. Wash step
 10. Air dry pellet.
 11. Resuspend pellet in the appropriate volume of TE or water once is fully dried.